Resorte SENTALLOY

SENTALLOY COIL SPRINGS

RESORTES SENTALLOY

- Open Coil Springs
Resortes abiertos
- Close Coil Springs
Resortes cerrados
- Sentalloy Coil Springs Comparative Studies
Estudios comparativos de resortes Sentalloy
- Differences Between Open and Close Coils
Diferencias entre resortes abiertos y cerrados

■ OPEN COIL SPRINGS

■ RESORTES ABIERTOS

Open coils springs are fabricated from Sentalloy wire, it possesses many desirable properties for physiologic tooth movement that are not possible with the stainless steel type coil springs, they deliver a gentle, continuous force to move teeth without periodontal stress, all that is needed is normal body heat, which activates the martensite transformation.

Open coils are 15 mm in length when passive and can be compressed to 3 mm. (Figure 3.41).

To engage the open coil, first mark the mesial and distal tooth destinations on the arch wire. No distortion will occur during dental movement, and there is limited patient discomfort (Figures 3.42 & 3.43).

Sometimes, when the open coil length is not enough, it is possible to use Sentalloy Stop Wound Open Coil Springs (Figure 3.44).

Those coils are 18 cm. in length, and sharp points are completely eliminated. After determining the

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
Fig. 3.42

Fig. 3.43

length of the coil needed, simply cut the wire at the center of the tight spiral as indicated (Figure 3.45).

Since the Sentalloy coil needs about 1/10 as much as steel, good hygiene is much easier, but do advise the patient to rinse and brush around the coil. No other adjustment is needed.

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CLINICAL APPLICATION

An open coil with 150 g. of superelastic activity is compressed between the maxillary central incisor and the canine. In three months the retraction of canine is completed (Figures 3.46 & 3.47).

Fig. 3.46

Maxillary molars can be moved distally in 4 to 7 months without patient cooperation using a modified Nance as an anchorage device. The coil is placed between the first bicuspid and first molar bypassing the second bicuspid. Cut the Open Coil so that they extend to the end of the buccal tube. Ligate the archwire to the first bicuspid. Next feed the arch wire into the tube, ensuring that the archwire is fully inserted with 2 mm. of wire extending beyond the tube. Uprighting springs may be used from the first bicuspid to distal of laterals to control any anterior movement (Figures 3.48 & 3.49).

Fig. 3.47

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
Fig. 3.48


Fig. 3.49

Figures 3.46 & 3.47 reprinted from Miura F, Mogi M, Ohura Y, Karibe M. The super-elastic property of the Japanese NiTi alloy wire for use in orthodontics. Part III Studies on the Japanese NiTi alloy coil springs. Am J Orthod Dentofac Orthop 1988;94:89-96 with permission from American Association of Orthodontists.

Figuras 3.46 y 3.47 reimpresas de Miura F, Mogi M, Ohura Y, Karibe M. The super-elastic property of the Japanese NiTi alloy wire for use in orthodontics. Part III Studies on the Japanese NiTi alloy coil springs. Am J Orthod Dentofac Orthop 1988;94:89-96 con permiso de la Asociación Americana de Ortodontistas.

■ CLOSE COIL SPRINGS

■ RESORTES CERRADOS

Sentalloy close coil springs deliver a gentle, continuous force to move teeth without periodontal stress. When warmed to body temperature the light continuous force and shape memory are activated.

Two designs are available: with Two eyelets or with an eyelet and an eyelet with ligature wire (Figure 3.50).

With two eyelets, the initial length is 3 mm. unactivated (This does not include the eyelet lengths). These can be activated to 15 mm without any deformation or change in force (Figure 3.51).

On an eyelet with ligature, you have the option of hooking the eyelet on the molar, then activating the spring and tying it to the cuspid bracket with the .008" ligature wire.

Should you wish to span an area greater than 15 mm., measure the span first, then place the coil on the molar hook and activate the coil 15 mm. finally tie the ligature wire to the bracket for the remainder of the span (Figure 3.52).

These eyelets are laser fused, which prohibits unwinding like other companies coils, which is why it is easier to use (Figure 3.53).

CLINICAL APPLICATIONS

Canine Retraction

A close coil spring with 120-130 g. of superelastic activity was attached for canine retraction,

Fig. 3.50

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Fig. 3.54

Fig. 3.55

and the distal ends connected to hooks on molars, in 5 months the distal canine retraction has been completed (Figures 3.54 & 3.55).

Anterior/Posterior Retraction

Springs are available for segmented retraction of the maxillary four incisors. The increased force in these springs was designed to overcome the resistance of the roots and to accomplish the same result as rectangular wire. You will experience better root parallelism control during the space-closing phase. Stops may be set on the arch wire to ensure the precise movement (Figure 3.56).

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Fig. 3.56

Figures 3.74 & 3.75 reprinted from Miura F, Mogi M, Ohura Y, Karibe M. The super-elastic property of the Japanese NiTi alloy wire for use in orthodontics. Part III Studies on the Japanese NiTi alloy coil springs. Am J Orthod Dentofac Orthop 1988;94:89-96 with permission from American Association of Orthodontists.

Figuras 3.74 y 3.75 reimpresas de Miura F, Mogi M, Ohura Y, Karibe M. The super-elastic property of the Japanese NiTi alloy wire for use in orthodontics. Part III Studies on the Japanese NiTi alloy coil springs. Am J Orthod Dentofac Orthop 1988;94:89-96 con permiso de la Asociación Americana de Ortodontistas.


Fig. 3.57. Giancoli A, Gianelly A. Three-dimensional control in extraction cases using a bidimensional approach. World J Orthod 2001;2:168-176.

1. 22 year old woman with bimaxillary protrusion and lip incompetency at rest.
2. Four first premolars were extracted.
3. 0.016 x 0.022-inch SS wire in maxillary arch with stops on maxillary first molars. Close Coil springs (150 g.) Class I forces from the molar hook to canines. Head gear was used for maximum anchorage and also used maxillomandibular Class II (100 g.) elastics from mandibular first molars to a hook on the maxillary arch.
4. Maxillary retraction was accomplished with sliding mechanics on a 0.018 x 0.022-inch SS. Mandibular incisors were retracted on a 0.016 x 0.022-inch SS wire. Close Coil springs exerting a force of 300 g. were extended from molar to distal of the lateral incisors.
5. Incisor retraction was completed in 3 months and the therapy lasted 24 months.

Fig. 3.57. Giancoli A, Gianelly A. Three-dimensional control in extraction cases using a bidimensional approach. World J Orthod 2001;2:168-176.

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SENTALLOY COIL SPRING COMPARATIVE STUDIES

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- Closed Coil Spring Comparative Studies
Estudio comparativo de resortes cerrados

- Open Coil Spring Comparative Studies
Estudio comparativo de resortes abiertos

- Differences between Open and Close Coils
Diferencias entre resortes abiertos y cerrados

Miura F, Mogi M, Ohura Y, Hamanaka H.

The super-elastic property of Japanese NiTi alloy use in orthodontics: Part III studies on the Japanese NiTi alloy coil springs. Am J Orthod Dentofac Orthop 1988;94:89-96.

CLOSED COIL SPRING COMPARATIVE STUDIES

A) Sentalloy springs and commercially springs

Comparative tensile studies of Closed Coil Sentalloy springs and the commercially available steel and Co-Cr-Ni coil springs were performed. The test coil spring size was 0.009 inch in diameter with a lumen of 0.030 inch. A tensile speed of 10 mm/min at a temperature of 37°±1°C was induced.

Findings

Load deflection curves of the commercial coils showed a linear relationship (Figure 3.58).

In contrast, Sentalloy closed coil springs showed an almost constant load and lower values than the first two samples. No permanent deformation occurred even when it was stretched to 500% of its original length. This is due to its unique superelastic property (Figure 3.59).

B) Diameter of the wire

The diameters of the wire selected for the testing were 0.009 inch, 0.010 inch and 0.012 inch. The lumen for all was 0.030 inch and martensite transformation temperature was at 0°C.

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


Findings

When the diameter of the wire increases while the lumen remains constant, the load value of the super-elastic activity increases but the range of super-elastic activity diminishes (Figure 3.60).

C) Size of the lumen

Coil springs tests were made from only 0.009 inch diameter with the lumen at 0.020, 0.027 and 0.030 inch and martensite transformation temperature was at 0°C.

Findings

When the size of the lumen increases while the diameter of the wire remains constant, the load value of the super-elastic activity decreases and the range of the super-elastic activity increases (Figure 3.61).

D) Martensite Transformation Temperature

Coil springs with three different martensite temperatures were examined: -43°C, -25°C and 0°C.

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

Findings

The load value of the superelastic activity is reduced with elevation of the martensite transformation temperature (Figure 3.62).

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

Melsen et al.²¹ provided data on load-deflection rate and maximal force generated for 19 coil springs. GAC close coil springs exhibited significant different activation and deactivation curves (Figure 3.63). Of the products tested only the springs from GAC exhibited the behaviour of the superelastic wire (Figure 3.64).

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
OPEN COIL SPRING COMPARATIVE STUDIES

A) Sentalloy springs and commercially springs

Comparative compression tests of stainless steel coil spring, Co-Cr-Ni coil spring and Sentalloy open coil spring samples were conducted on the 50 mm. open coil spring. A core was placing in the coil and its size was one half the size of the lumen. The test coil spring size was 0.009 inch in diameter with a lumen of 0.030 inch. A compression speed of 10 mm/min at a temperature of 37°±1°C was induced.

Findings

The results of the compression test showed that commercially available open coil springs demonstrated a linear relationship illustrated by the load deflection curves. After the coil loops were compressed tightly into tight contact with each other, permanent deformation occurred (Figure 3.65).

In contrast, Sentalloy open coil springs demonstrated a deflection curve with a constant load value and no permanent deformation occurred, and a light and continuous force value was exhibited (Figure 3.66).

B) Diameter of the wire

The diameters of the wire selected for the testing were 0.009 inch, 0.010 inch and 0.012 inch. The lumen for all was 0.030 inch and martensite transformation temperature was at 0°C.

Findings

When the open coil is compress, the load value of the super-elastic activity increases in proportion to the increased size of the diameter of the wire (Figure 3.67).

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

C) Size of the lumen

Coil springs were made from only 0.009 inch diameter with the lumen at 0.020, 0.027 and 0.030 inch and martensite transformation temperature was at 0°C.

Findings

When the size of the lumen increases while the diameter of the wire remains constant, the load value of the super-elastic activity decreases and the range of the super-elastic activity increases (Figure 3.68).

D) Martensite Transformation Temperature



Coil springs with three different martensite transformation temperatures were examined: -34°C, -15°C and 0°C.

Findings

When the martensite transformation temperature is elevated the load value of super-elastic activity is reduced. It is possible to change the load value of super-elastic activity in open coil springs having the same sizes of wire diameter and lumen. Between these points the coils will exert 100 (light), 150 (medium), or 200 (heavy) grams of force until 14 mm, and then the force falls to 0 as it becomes passive again (Figure 3.69).

Conclusion

Studies has shown that the load value of superelastic activity of Sentalloy Open and Close coil springs can be controlling changing the specifications of the wire such as the diameter

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lumen, pitch of coils and martensite transformation temperature. It is now possible to fabricate coils, springs that possess specific load values and ranges of super-elastic activity for specific tooth movement.

DIFERENCES BETWEEN OPEN AND CLOSE COILS

Both of them are made from Sentalloy wire, but due to some differences during activation and deactivation phases they demonstrate some important differences.

With a 0.009 inch diameter and a 0.030-inch sample, open coil springs show a more constant load value of superelastic activity than the closed coil springs. The closed coil spring shows a value of 100-200 g. (Figure 3.70) but the open coil exhibits a load value of 100g. (Figure 3.71).

This is because during the deactivation phase, when the close coils are stretched, the lumen size becomes smaller and the inclination of load deflection curve becomes acute (Figure 3.70).

On the other hand in the open coil springs when they compressed, the lumen size becomes larger and the inclination of the curve will be obtuse (Figure 3.71).

That is why a more desirable continuous force can be obtained from the open than the closed coil spring.

Terminal period load deflection curve

Because the super-elastic activity decreases to 0 in the final stage, sometimes the load value is not enough for tooth movement, in case of open coil it

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


3

should be a little longer than the desired distance and the coil should be compressed tightly at the time of application of the pressure. In case of closed coil the activated length should be a little longer than the required distance and the tension should be just before the permanent deformation stage. By means of fabricating the new closed coil, the individual coil-to-coil interface area, it possesses a definite force activity that compel the coils to tend to remain constantly in tight contact with each other, thus tooth movement continues until the very end of the deactivation stage (Figure 3.72).

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Open Coils, Resortes Abiertos

Giannelly A, Bednar JA, Dietz VS. Japanese NiTi coils used to move molars distally. *Am J Orthod* 1991;99:564-566.

Karibe M. Orthodontic studies on super-elastic NiTi alloy wire coil spring. *J Jpn Orthod Soc* 1991;50:61-74.

Samuels RH, Rudge SJ. A quick and simple method of reactivating space-opening push-coil springs. *J Clin Orthod* 2000 December;XXXIV(12):706-708.

Closed Coils, Resortes Cerrados

Andres LS. Comparison of NiTi coil springs vs. elastics in canine retraction. *J Clin Orthod* 1994 May;29:293-295.

Küseler A, Melsen B. Orthodontic treatment: art or science? *J Clin Orthod*;34(8):371-375.

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Samuels RH, Rudge SJ, Mair LH. A comparison of the rate of space closure using a nickel-titanium springs and an elastic module: a clinical study. *Am J Orthod Dentofac Orthop* 1993;103:464-467.

Samuels RH, Rudge SJ, Mair LH. A clinical study of space closure with nickel-titanium closed coil springs and an elastic module. *Am J Orthod Dentofac Orthop* 1998;114:73-79.

Figures 3.58-3.62 & 3.65-3.72 reprinted from Miura F, Mogi M, Ohura Y, Karibe M. The super-elastic property of the Japanese NiTi alloy wire for use in orthodontics. Part III Studies on the Japanese NiTi alloy coil springs. *Am J Orthod Dentofac Orthop* 1988;94:89-96 with permission from American Association of Orthodontists.

Figuras 3.58-3.62 y 3.65-3.72 reimpresas de Miura F, Mogi M, Ohura Y, Karibe M. The super-elastic property of the Japanese NiTi alloy wire for use in orthodontics. Part III Studies on the Japanese NiTi alloy coil springs. *Am J Orthod Dentofac Orthop* 1988;94:89-96 con permiso de la Asociación Americana de Ortodontistas.